


ENDO ACTIVATOR[®]

Technique Card



PRODUCT APPLICATIONS

- Debridement and disruption of the smear layer and biofilm
- Placement of calcium hydroxide and MTA around root curvatures
- Removal of residual obturation materials during retreatment

ACTIVATE YOUR ENDODONTICS

PRODUCT FEATURES

EndoActivator Tips

- Easy snap-on/snap-off design
- Strong, flexible, medical-grade polymer composition
- Color-coded by size
- 22 mm lengths with depth gauge rings at 18, 19 and 20 mm

Handpiece / Driver

- Contra-angled design for easy access to posterior teeth
- 3-speed sonic motor for various clinical applications.
- Ergonomic, cordless, and battery-operated handpiece



U.S. Patent Nos. 7,261,561, 8,235,719, 8,328,552 and 8,388,345
International Patent Nos. 1,743,593 and 1,712,201
Additional Patents Pending

ENDO ACTIVATOR[®]

SAFE • EFFECTIVE • SIMPLE

TECHNIQUE SEQUENCE*

- Prepare canal to produce a fully tapered shape
- Fill pulp chamber with EDTA, QMix or other final rinse solution
- Select the Activator that manually fits loosely within 2 mm of working length
- Place the barrier sleeve over the Driver to protect the entire handpiece
- Attach the Activator over the barrier-protected Driver
- Place the attached Activator tip into the prepared root canal
- Depress the ON/OFF switch to activate solution
(Note: Switch defaults to 10,000 cpm, as indicated for cleaning procedures)
- Use a pumping action to move the Activator in short, 2-3 mm vertical strokes
- Hydrodynamically agitate the intracanal solution for 60 seconds
- Irrigate with NaOCl, then use intracanal suction to eliminate loose debris

**Sequence for cleaning root canal systems. For other applications, see Directions For Use.
For scientific evidence supporting the use of the EndoActivator, visit endoruddle.com*



See Directions For Use

ADVANCED  ENDODONTICS[®]

Clifford J. Ruddle DDS FACD FICD

(800) 753-3636 • (805) 964-8838 • www.endoruddle.com