

ACTIVATE YOUR ENDODONTICS

The EndoActivator™ is designed to safely and vigorously energize the hydrodynamic phenomenon. Evidence-based endodontics has shown that cavitation and acoustic streaming significantly improve debridement and the disruption of the smear layer and biofilm. Activated fluids promote deep cleaning and disinfection into lateral canals, fins, webs, and anastomoses. A cleaned root canal system facilitates 3-D obturation and long-term success.

Activate your endodontic solutions today with this easy-to-use technology!




DESIGNED BY
Clifford J. Ruddle DDS
Robert H. Sharp DDS
Pierre Machtou DDS



EndoActivator Tips
(Shown Actual Size)

CLINICAL FEATURES

- Ergonomic handpiece design
- Strong, flexible polymer tips
- Effective deep lateral cleaning
- Safe and efficient technology
- Simple, affordable system



EndoActivator Handpiece/Driver
(Shown Actual Size)

(800) 753-3636

www.endoruddle.com

U.S. and International Patents Pending


ENDO ACTIVATOR™

SAFE • EFFECTIVE • SIMPLE


SYSTEM SIMPLICITY

The EndoActivator Tips are 22 mm in length and are available in three sizes: small, medium, and large. Specific features include:

- Easy snap-on/snap-off design
- Strong, flexible, medical-grade polymer composition
- Clean guard component for the deflection of aerosols
- Color-coded by size for easy identification
- Convenient depth gauge rings at 18, 19 and 20 mm
- Sold separately in reorder packs or available as a complete system with the Handpiece/Driver



Handpiece/Driver



The EndoActivator System is designed to provide a safer, better, and faster way to achieve success in the following clinical applications:

- Debridement and disruption of the smear layer and biofilm
- Placement of calcium hydroxide and MTA around root curvatures
- Removal of residual obturation materials during retreatment procedures

ADVANCED  ENDODONTICS®

Clifford J. Ruddle DDS FACD FICD
227 Las Alturas Road Santa Barbara, CA 93103
(800) 753-3636 • (805) 964-8838 • Fax (805) 965-8253

(800) 753-3636
www.endoruddle.com

U.S. and International Patents Pending

